Upcoming entrained bubbles get eliminated through a holding time of approx. 1 hour with FotoDent® model.

To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and colour variations of the FotoDent® model material.

Possible dangers

Please note the safety instructions for processing FotoDent® model.

Hazard statements: H317 May cause an allergic skin reaction. H412 Harmful to aquatic life with long lasting effects.

Precautionary statements: P261 Avoid breathing dust/fume/gas/mist/vapours/spray. P273 Avoid release to the environment. P280 Wear protective gloves/protective clothing/eye protection/face protection. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P333+P313 If skin irritation or rash occurs: Get medical advice/attention. P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains: 1,4 Butanediol dimethacrylate, Diurethandimethacrylate, hydroxypropylmethacrylate, monoester with propane-1,2-diol, Diphenyl(2,4,6-trimethylbenzoyl)phosphinoxide

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® model packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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Product description

FotoDent® model is a light-curing resin for the production of dental models using 385 nm / 405 nm LED based stereolithography. For the post curing of FotoDent® model built parts, the light curing units FotoDent® flash and PCU LED are recommended, which are optimal balanced to the production process.

Processing

- During the application of FotoDent® model we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- FotoDent® model should be well shaken approx. 1 hour before use.
- Carefully pour FotoDent® model into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.
- Processing temperature:
Rapid Shape printer: $23 \pm 3^\circ\text{C}$
Asiga Freeform Pro2 UV: $35 \pm 3^\circ\text{C}$
- Select the product specific parameter settings for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 10 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended.
- Post-curing:
FotoDent® flash 12 minutes + 1 minute gas pre-flow
PCU LED 8 minutes, 80 % light power (vacuum)
- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® model. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content.