Dreve Design Rev.160118-DU190118-500-39342

Upcoming entrained bubbles get eliminated through a holding time of approx. 1 hour with FotoDent® tray.

To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and colour variations of the FotoDent® tray material.

Possible dangers

Please note the safety instructions for processing FotoDent® tray.

Hazard statements: H318 Causes serious eye damage. H317 May cause an allergic skin reaction. H412 Harmful to aquatic life with long lasting effects.

Precautionary statements: P261 Avoid breathing dust/fume/gas/mist/vapours/spray. P273 Avoid release to the environment. P280 Wear protective gloves/eye protection/face protection. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P310 Immediately call a POISON CENTER/doctor. P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains: 2-Propenoic acid, reaction product with Pentaerythrite, 1,4 Butanediol dimethacrylate, phenyl bis(2,4,6-trimethylbenzoyl)-phosphineoxide, Pigments

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® tray packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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Intended use

FotoDent® tray is a light-curing resin for the production of dental impression trays using 385 nm LED based stereolithography. For the post curing of FotoDent® tray built parts, the light curing units FotoDent® flash and PCU LED N₂ are recommended, which are optimal balanced to the production process.

Processing

- During the application of FotoDent® tray we recommend wearing personal protective equipment e.g. suitable gloves, safety goggles etc.
- FotoDent® tray should be well shaken approx. 1 hour before use.
- Carefully pour FotoDent® tray into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.
- Processing temperature:
Carbon M1 / M2: 23 ± 3 °C
- Select the product specific parameter settings for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 10 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended.
- Post-curing:
PCU LED N₂ (vacuum) 8 minutes, 80% light power
- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® tray. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content.