

## 8. Hazards identification

Hazard and safety information can be found on the product label and/or in the corresponding safety data sheet.

## 9. Disposal

Disposal of contents/container in accordance with local/regional/national/international regulations and according to the safety data sheet.

## Glossary

|  |                                  |   |  |
|--|----------------------------------|---|--|
|  | Use-by date                      | Rx only   | Use only by professionals  |
|  | Please note instructions for use |  | Protect from sunlight  |
|  | Temperature limitation           | MD  | Medical device   |
|  | Manufacturer                     | CE  | CE marking of conformity according to the requirements for medical devices |
| LOT  | Batch code                       |   |  |
| REF  | Item number                      |  | Made in Germany  |

## Form of delivery

REF D353002

FotoDent® tray2 green transparent

1.0 kg bottle



Rx only

MD



# FotoDent® tray2

Date of issue: 2023-06-22/Rev. 00

## Directions for use



Dreve Dentamid GmbH  
Max-Planck-Str. 31 • 59423 Unna/Germany  
www.dentamid.dreve.de



CE



## Intended use

Light curing material for production of dental impression tray.

**Indications for use:** Manufacturing dental impression trays. The dental impression tray will be filled with dental impression material to reproduce the structure of a patient's teeth and gums.

**Indication:** Individual impression tray for taking dental impressions.

**Contraindication:** If an acrylic allergy exists, an allergic reaction may occur.

**Patient target group:** Patients who need dental care.

**Intended users:** The dental technician or specialist staff in the dental laboratory is responsible for the production of the impression tray in 3D printing. The dentist uses the finished impression tray on the patient.

## 1. Requirements / equipment

**Printer:** Carbon M2, M3, M3 Max

**Software:** Carbon Print Software

**Cleaning:** Isopropanol, orbital shaker, container

**Post curing unit:** PCU LED N<sub>2</sub>

## Basic Material Properties

Please refer to Technical Data Sheet (TDS). Available on request.

## 2. Process description

### 2.1 Preparation

- During processing of the liquid material, we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- Carefully pour the material into the pre-determined container of the production unit.
- Please see further steps within printer manual.
- Remove all bubbles with a cleaned object.

### 2.2 Printing Parameters

- Layer thickness: 300 µm
- Minimum wall thickness impression tray: 2.0 mm
- Minimum wall thickness impression handle: 5.5 mm
- Optimal Orientation: the handle faces the building platform and is built directly onto the building platform
- Supports: optional, to be added by the user

### 2.3 Printing process

- Processing temperature: 22 ± 3 °C
- Select the corresponding profile for FotoDent® tray2 in Carbon's software. Make sure the software used is up to date.
- Start the printing process.

### 2.4 Processing after the printing process

- It is recommended that the printing result be further processed directly after printing (see following steps).
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.

### 2.5 Cleaning

- Clean in isopropanol medium (purity ≥ 97 %) at an orbital shaker. Fill two suitable containers with isopropanol.
  1. Pre-clean with isopropanol for 6 minutes in the first suitable container.
  2. Blow off using compressed air.
  3. Perform the main cleaning step with isopropanol for 6 minutes in the second suitable container.
  4. Blow off using compressed air.

### 2.6 Post-curing

PCU LED N<sub>2</sub>      8 minutes, 80 % light output using nitrogen

## 3. Important notes

- The specifications and biological safety have been qualified using the aforementioned printer, the associated software and the process parameters indicated.
- To avoid deterioration of the material quality, protect the

liquid material from exposure to light. Variations from the described manufacturing process may impair biological safety, lead to modified mechanical characteristics and/or color variations of the material.

- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. However, thanks to the low viscosity, it is possible to filter FotoDent® tray2. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hours with FotoDent® tray2.

## 4. Cleaning and disinfection

For cleaning and disinfection, we recommend the use of a 0.2 % chlorhexidine solution or 0.5 % sodium hypochlorite solution with an exposure time of 5 minutes for disinfection, and then allow to dry completely.

## 5. Sterilisation

FotoDent® tray2 cannot be sterilised.

## 6. Serious incidents

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

## 7. Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® tray2 packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date