PCU 90

A Dreve Company

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EN Instructions for use

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Contents (EN)

1	Decommissioning	20
2	Unit description	20
3	Technical data	20
4	Declaration of conformity	21
5	Safety notes	22
6	Start up	24
7	Information on using the shielding gas	24
8	Functional elements	25
9	Operating field and control	25
10	Operation	27
10.1	Post curing at the upper level	27
10.2	Post curing at the lower level	28
10.3	Rotary plate post curing	28
10.4	Display and resetting of the operating hours	28
11	Bulb change	28
12	Care and maintenance	29
13	General notes	29
14	Fault rectification plan	30
15	List of spare parts	32
16	Wiring diagram	48–53

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deutsch
english
français

1 Decommissioning

The DIRECTIVE 2012/19/EU OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 4 July 2012 on waste electrical and electronic equipment (WEEE) lays down measures to protect the environment and human health. Waste electrical and electronic equipment marked as such may not be disposed of as unsorted municipal waste within the European Union. For detailed information regarding the professional disposal of WEEE please contact Dreve or your distributor.

2 Unit description

The PCU 90 is an industrial light polymerisation unit for curing high-quality components manufactured using stereolithography, DLP (digital light processing) or LCD (liquid crystal display) processes.

Thanks to the metal halide lights as well as the tuning of the wavelength, the PCU 90 is ideally suited for the post-curing of light-polymerising resins.

To avoid the dispersion layer on the surface of the printed components, the PCU 90 enables curing under the usage of inert gas.

The extremely large polymerisation chamber enables the curing of many large components at the same time

The control of the PCU 90 and the intuitive menu management makes operation of the unit straightforward. Only two parameters are set: the curing time, on a large operating button and the gas flow rate on the small setting button with locknut. The air current was adapted to requirements with baffles in order to optimise the cooling efficiency. A temperature sensor protects the unit and it will switch off the post curing process following overheating in an emergency.

3 Technical data

H x W x D	725 mm x 1185 mm x 612 mm
Polymerization chamber	360 mm x 610 mm x 410 mm
Weight	140 kg
Mains voltage	230 V / 50 Hz, 220 V / 60 Hz, CEE (Voltage / frequency)
Power consumption	3500 W
Protection	16 amps anti-surge
Maximum operating pressure	5 bar (70 psi)
Gas connection	Plug and socket connection for 8 mm x 6 mm x 1 mm hose

4 Declaration of conformity

Low Voltage Directive (2014/35/EU)

Herewith we declare that the machine described below in its conception and design and in the shape delivered by us is in accordance with the fundamental requirements for safety and health as prescribed in the applicable EC directives.

In the case of any change or modification of the machine, not authorised by us, this declaration becomes invalid.

Denomination of the machine: PCU 90
M43000 / M43000B
Unit type: Light polymerization unit

Applicable EC directives:
2014 / 30 / EU EMV Directive
2006 / 42 / EC According to Machinery Directive
2014 / 35 / EU Low Voltage Directive

Applicable EC directives:

2014 / 30 / EU
EN 61000-3-2:2014

EMV Directive

Electromagnetic compatibility (EMC) – Part 3-2: Limits – Limits for harmonic current emissions (equipment input current ≤ 16 A per phase) (IEC 61000-3-2:2014); German version EN 61000-3-2:2014.

EN 61000-3-3:2013

Electromagnetic compatibility (EMC) – Part 3-3: Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current ≤ 16 A per phase and not subject to conditional connection (IEC 61000-3-3:2013); German version EN 61000-3-3:2013.

EN 61326-1:2013

Electrical equipment for measurement, control and laboratory use – EMC requirements – Part 1: General requirements (IEC 61326-1:2013); German version EN 61326-1:2013. Emission according to living area, business and industrial undertakings as well as small enterprises. Interference stability according to industrial area.

2006 / 42 / EC

EN ISO 12100:2010
+ Corrigendum to
DIN EN ISO 12100:2011-03

Machinery Directive

Safety of machinery – General principles for design – Risk assessment and risk reduction (ISO 12100:2010); German Version EN ISO 12100:2010 +Corrigendum to DIN EN ISO 12100:2011-03.

EN 60204-1:2006 / A1:2009
+ Corrigendum to
EN 60204-1:2006

Safety of machinery – Electrical equipment of machines – Part 1: General requirements (IEC 60204-1:2005 / A1:2008); German Version EN 60204-1:2006 / A1:2009 + Corrigendum to EN 60204-1:2006; German Version CENELEC-Cor.:2010 to EN 60204-1:2006.

2014 / 35 / EU

EN 61010-1:2010

Low Voltage Directive

Safety requirements for electrical equipment for measurement, control and laboratory use – Part 1: General requirements (IEC 61010-1:2010 + Cor.:2011); German version EN 61010-1:2010.

Innovation MediTech GmbH
Max-Planck-Str. 31 · 59423 Unna/Germany
Tel.: +49 2303 8807-0 · www.innovation-meditech.de

5 Safety notes

Caution! Read these notes through carefully before connecting and starting up the unit. The operational safety and the functioning of the unit can only be guaranteed if the general regulations on safety and prevention of accidents of the legislator and the safety notes in the operating instructions are observed.

1. The unit may only be used in line with these operating instructions. We will not be held liable for any damage which may result from inappropriate use or incorrect operation.

2. The voltage stated on the rating plate must match the voltage of the power source.

3. Only operate the unit from a socket with a protective earthing conductor. Never touch the power plug with wet hands.

4. The unit must stand on a level and firm foundation which is sufficiently stable and durable for the 140 kg approximate operating weight.

5. Do not introduce any objects into the unit from outside the polymerization chamber.

6. Any accessories and spare parts used must be cleared by the manufacturer. We will not accept any liability for damage incurred due to the use of non-standard parts.

7. Do not lock the operating and safety elements of the unit in place.

8. The good condition and operational safety of the unit must be checked prior to each use. If the condition is not fault free the unit must not be used and it must be appropriately marked.

9. The power plug must be pulled out prior to any cleaning and maintenance work on the unit or the exchange of parts.

10. Signs and stickers must be kept in legible condition and they must not be removed.

11. Only authorized specialists may open the unit up and carry out maintenance work.

12. Only nitrogen or carbon dioxide may be used as protective gases. Do not use any flammable or poisonous gases. Always observe the general regulations for safety and prevention of accidents when handling the chosen gas.

13. Only use approved gas bottles.

14. Check the supply hose for leak tightness when connecting the gas bottle to the unit. Make sure that the gas bottle cannot tip over.

15. Never connect the gas bottle without using a pressure regulator. The unit may be operated with a maximum pressure of 5 bar (70 psi).

16. Shut off the gas supply if the unit is not being used. Never operate the unit unsupervised.

17. Ensure that the work area has adequate ventilation.

18. The pressure of the gas line must not exceed 5 bar (70 psi).

19. Do not keep any flammable substances close to the unit.

20. Approved user: the company which operates the unit must make the operating instructions available to the user and ensure that he has read and understood them. The user may not operate the unit until this has been done.

21. In order to prevent water entering the unit (e.g. spray) should be set up in a dry area.

22. The unit must be disconnected from the mains if it is not being used.

23. Only touch the front flap by the handles in order to prevent squashing.

24. Do not operate the unit if the assembly flap on the underside of the unit is open.

25. Unauthorized conversions and changes are not permitted for safety reasons.

26. The operating and maintenance conditions in these operating instructions must be strictly adhered to. When working with the PCU 90, the general regulations for prevention of accidents must be observed.

27. Do not touch the glass slabs and rotary plates as they become very hot during operation.

28. Do not obstruct or block the ventilation slot on the left hand side of the unit. Always leave sufficient distance to the wall.

6 Start up

Unpack the unit and place it on a sufficiently stable substructure (load bearing capacity at least 150 kg). Remove film carefully and examine the unit for signs of transport damage. Make CEE electrical connection. Open flap to curing area. Pull glass insert forwards with metal frame until the package can be removed. Remove gas hose from the package and check the rotary plate for damage. Remove the three upper Styrofoam inserts. Take the drawer right out. Remove the protective film from the displacer. Connect

the gas hose and set the pressure regulator on the gas bottle to 5 bar.

Supply scope

- 1x PCU 90
- 1x 8 mm gas hose line incl. quick coupling
- 1x glass insert
- 1x rotary plate made of special glass
- 1x operating instructions
- 1x displacer (displacement body for volume reduction, reduced gas consumption)


7 Information on using the protective gas

Only nitrogen or carbon dioxide may be used as protective gases. Do not use any flammable or poisonous gases. Always observe the general regulations for safety and prevention of accidents when handling the chosen gas. When using the protective gas it is possible to carry out the polymerization in an atmosphere with little oxygen. The protective gas replaces the air in the polymerization chamber and thus reduces the oxygen content in the polymerization area of the PCU 90.


This means that it is possible to post cure components without an inhibition layer. The protective gas supply is controlled automatically once the polymerization process starts.


The flow rate is pre-set (30 l/m). The adjustment can be changed on the regulating needle valve on the front of the unit (observe the safety notes). The output pressure of the bottle can be modified on the pressure regulator of the gas bottle and it must be set to 5 bar.

8 Functional elements

9 Operating field and control

10 Operation

Secure CEE power supply. Switch on the main switch low down on the right of the operating side of the unit. Connect protective gas. Set gas bottle pressure to 5 bar. The following will appear on the display:

0:00

The lighting time is set by turning the multifunction knob. From 10:0 the seconds are displayed in stages of 30.

10:0 10:3 11:0

It is not possible to set more than 20:0 of lighting time. Pressing the button starts the process. A long press (3 seconds) on the button always stops the respective program.

The PCU 90 can work with two different lighting concepts: the area lighting with six lamps which lights up the lighting field optimally from above and also the rotation function which provides optimum lighting of the area above the rotary plate with four lamps from above and two from the side.

The setting of the gas flow is performed at the setting button on the right. Release the locknut, set the desired gas flow and retighten the locknut. The gas flow is set at approx. 30 l / m in the factory for post curing at the higher level.

10.1 Post curing at the upper level

Move the flip switch upwards to area. Open flap to curing area. Remove the rotary plate from the lower level if it is present. Remove the insert and put the displacer in position.

Caution: Never let go of the opening flap. It must always be completely guided by your hand whilst opening and closing. Push the glass slab in. Put parts for post curing into position and close opening flap. The magnetic locking is audible. Set post curing time. A quick press of the multifunction button starts the post curing program and the set post curing time flashes. Set the gas

flow to approx. 30 l/min with the help of the setting screw. The value is shown on the blue flow display. A program for pre-gassing and filling the lighting chamber and protecting the lights will run for four minutes. After this the actual lighting period starts. The display expires. After three minutes the gas valve closes automatically and the blue flow display goes out so that the PCU 90 is operated with protective gas for seven minutes.

Normally an oxygen content is achieved in the chamber which is needed for inhibition free post curing. Once the

Program has finished a short bleep will be audible and the display will return to the previously set value. The components can then be removed. Functional sequence as for the upper level. Please remove the glass insert

and the displacer as well, modify the gas flow to the greater volume and set to at least 50 l/m. Push the glass insert in at the lower level, cover with parts for post curing and then continue as described in 10.1.

10.3 Rotary plate post curing

Move the flip switch upwards to rotation. Open flap to curing area. If present, remove the insert and the displacer. Push the rotary plate in carefully. Put parts for post curing into position and close opening flap. The magnetic locking is audible. Set post curing time. A quick press of the multifunction button starts the post curing program and the set post curing time flashes. Set the gas flow to approx. 50 l/min with the help of the setting screw. The value is shown on the blue flow display. A program for pre-gassing and filling the lighting chamber and protecting

the lights will run for four minutes. After this the actual lighting period starts. The display expires. After three minutes the gas valve closes automatically and the blue flow display goes out so that the PCU 90 is operated with protective gas for seven minutes.

Normally an oxygen content is achieved in the chamber which is needed for inhibition free post curing. Once the program has finished a short bleep will be audible and the display will return to the previously set value. The components can then be removed.

10.4 Display and resetting of the operating hours

In order to guarantee better control of the output and service life of the bulbs, the PCU is equipped with an operating hours counter. If the multifunction button is pressed and the machine's main switch is activated, the operating hours

will be shown. After that keep the button pressed down in order to return to the lighting menu. The operating hours counter can only be reset by an authorized service technician.

11 Bulb change

A bulb change should only be carried out by an authorized service technician.

The instructions for changing bulbs are with the exchange lights.

12 Care and maintenance

Maintenance work to be carried out by the user

Maintenance unit	Frequency of maintenance
Clean glass screens of the unit with isopropyl	daily
Before start up: Check whether all connections of the unit are seated correctly: <ul style="list-style-type: none"> • Gas connection • Power connection 	daily
Cleaning of the housing with a mild cleaning agent	weekly
Check whether all seals are functioning properly and/or are seated correctly	monthly
Bulb change once the service life of the lights has reached > 750 h.	once the operating hours counter shows it has been reached (shown on display)
The braked flap holders must regularly be checked for damage and wear. When adjusting the braking force, make sure that the adjusting screws of both flap holders are evenly adjusted.	monthly

13 General notes

Repair and maintenance work on the PCU 90 may only be undertaken by specialized personnel. The power feed to the unit must be disconnected in each case. Innovation MediTech GmbH will only be responsible for the safety, reliability and performance of the unit if

- the unit is only used for the appropriate purpose.
- it is only operated by specially trained staff.

- the unit is used in accordance with these operating instructions.
 - modifications, new settings, changes or repairs are only carried out by personnel that have authorized it.
- As we are constantly working on further development of our products, we reserve the right to make technical changes.

14 Fault rectification plan

Fault	Cause	Remedy	Fault	Cause	Remedy	
Appliance does not show any function	Power connection incorrect	Check power connection	Parts are sticky	Gas bottle empty	Exchange gas bottle	
	Main switch off	Switch main switch on		Pressure regulator is closed	Open pressure regulator	
	Safety device has triggered	Switch safety device on		Flap seals defective	Renew flap seals	
Main switch on, time circuit does not show any function	Main switch defective	Exchange main switch		Operating pressure on pressure regulator is no longer at 5 bar	Set pressure regulator correctly to 5 bar	
	Power board defective	Exchange power board		Regulating needle valve was set differently	Set regulating needle valve correctly	
Time circuit functioning	Connection cable from the power board to control board has been released	Connection cable must be inserted properly		Regulating needle valve defective	Exchange regulating needle valve	
	Flap not closed	Close flap		Relay board defective	Exchange relay board	
Program cannot be started	Flap switch defective	Exchange flap switch		Magnetic valve defective	Exchange magnetic valve	
Objects become too warm	Ventilator fan defective	Exchange ventilator fan		Operating flap is not locked	Program must be started beforehand	Start program
					Flap is not closed centrally and properly	Close correctly with both handles
	Electromagnets are defective	Exchange electro-magnets				
	Relay board defective	Exchange relay board				

15 List of spare parts:

	REF
1 PCU M120 housing	540400
2 MH Module 250 W XL	80034746
3 UV Lamp – Lamp excluding Reflector	5010190
4 Reflector MH 250 Module	5010151
5 Ballast MI 530-250 TS	540438
6 Ignitor ELI400	540439
7 Mini switching power supply, enclosed 24 V/1.1A 25 W	540405
8 Borosilicate screen 605 x 405 x 3mm	540407
9 Borosilicate screen D = 350 x 3mm	540408
10 Borosilicate screen 625 x 430 x 3mm	540409
11 Borosilicate screen 280 x 170 x 3mm	540410
12 Brake flap stay left with fine setting	540412
13 Brake flap stay right with fine setting	540413
14 Disc magnet	540416
15 Flow meter SFAB SA-M12 2:4	540418
16 Grommet HV 1302.	50892
17 Normplast ZEV 6–9/PG9 white	50869
18 Drufosmart electronics	55201X4
19 Chrome plated rotary knob D=32mm	51322
20 Keyboard film Drufosmart basic	55820
21 Housing foot 30/25-SW	51166
22 Axial ventilation fan 4650N,	50909
23 Werit Porcelain terminal block 1204-2	51096
24 Hinge 50 x 50 mm	50651
25 Distance sleeve 19.5 x 15.5 mm	55978
26 Microswitch V-153-1C5	50963
27 Magnetic clamp HT-D 25-F 24VDC 100 %ED	5755
28 Cover handle for Polyflux 2000	50947
29 Isofix P-profile seal white	51206
30 Cellular rubber-square line	57027
31 Small contactor relay 4S UC230 V 50 Hz	608052
32 Small contactor relay E3250, 3250-00310	60805
33 Cimco mounting pad 27 x 27 mm	5588
34 Werit 1204 Porcelain 3-pin.	51158
35 Werit elastic plug terminal	51063
36 Insulation distance sleeve, 3/6 x 7.5 PA	50959
37 Pur plastic hose type: H-PU-8/6-blue	50341
38 Norma chuck jaw clamp steel galvanised 1part	2006

REF

39 Restrictor valve, # 6.01.18/1N, G1/8	50937
40 3/2 way magnetic valve – 24V DC	52021-24
41 Angled screw in connection mini	5533
42 Hose, blue	5537
43 Push in fitting type: B-GE-8-1/8 MSv-bl	50340
44 Connection cable	540419
45 3-way distributor plug connection 6/4 G1/4"	540420
46 Screw in plug connection 90°	540421
47 Screw in connection G1/4" 8/6 hose	540422
48 Plug in connector 8 hose both sides	540423
49 Screw in connection G1/8" 8/6 hose	540424
50 Screw in plug connection 90° G1/8" 8/6	540425
51 Plastic rotary plate holder	540428
52 Drive motor	540429
53 Mains filter Schaffener FN2090-16-06	540432
54 Brass sleeve QM 1/8-1/8	50269
55 Anti-interference capacitor 0.047yF, 275 V ac 11 x 18 x 5.5D	50799
56 CE connector blue	
57 Thermostat	50759
58 Red warning light	50857
59 Time relay	57029
60 Safety device	
61 Stainless steel clamp	51201
62 Pad for polymerization chamber	M43013
Only for M43000B 220V / 60Hz	
3 Ballast MI 620-250 TS	540442
29 Small contactor relay 4S UC220 V / 60 Hz	60805260
30 Small contactor relay BR010 Uc 220 V / 60 Hz	540443

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La société Innovation MediTech GmbH n'est responsable de la sécurité, fiabilité et performance de l'appareil que si