Bedienungsanleitung
Instructions for use
Mode d'emploi

deutsch

english

français

PCU LED

Pioneering Technology

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**Light polymerization device for the
CAD/CAM 3D printing process**

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Contents (EN)

1. Device description	page	24
2. Technical data	page	25
3. Declaration of conformity	page	26
4. Safety instructions	page	27
5. Commissioning	page	29
6. Information on the use of the protective gas	page	29
7. Functional elements	page	30
8. Operation	page	31
9. Maintenance	page	37
10. Monitoring	page	38
11. Troubleshooting plan	page	38
12. Exploded view	page	40
13. Spare parts list	page	41
14. General information	page	42
15. Warranty	page	43
16. Decommissioning of the product	page	43

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1. Device description

The PCU LED (Post Curing Unit LED) is a light polymerization device for the curing of components which have been created by means of 3D printing. The following device properties are used in the curing process:

- a) LED technology
- b) Electronic control (micro-controller)
- c) Logging and monitoring of the process parameters
- d) Simple control concept
- e) Process reliability
- f) Different curing environments (air, vacuum, nitrogen)
- g) Open system

Device properties:

a) LED technology
The advantages of the LED technology are that a permanent light of one wavelength is emitted. The permanently emitted light radiation can be captured by a standard UV sensor and thus allows for the measurement and the observation of the light output in the exposure chamber. Furthermore, LED lamps do not produce any heat which could cause the material properties of the plastic component to be cured to be negatively influenced. Another advantage of the LED technology is that the stated operating life of the LED lamps is approx. 20,000 hours and is accompanied by low operating and maintenance costs.

b) Electronic control (micro-controller)
The use of a micro-controller allows for

the control and reading of the current, pressure and temperature sensors.

c) Logging and Monitoring

The information which is recorded by the current, pressure and temperature sensors represents the used process parameters (e.g. light on/off; vacuum on/off; nitrogen on/off). These process parameters are recorded, saved and transferred to a USB storage medium. This allows for the transparent traceability of the production chain of the medical device.

d) Simple control concept

The device is operated by the removable rotary button. The corresponding option is selected by turning the button. The selected option is operated by pressing. Longer pressing opens the parent and sub-folder menu. Another removable element is the handle which is used to open the lid of the exposure chamber. Both elements can be cleaned with the appropriate cleaning agents after removal.

e) Process safety

To be able to guarantee a high process liability, the control software is subdivided into two levels. Level 1 is accessible to the user and process owner. Level 2 is only accessible to the process owner. In level 1 only a predefined process can be selected. In level 2, for example, a process can be defined and an option can be switched on and/or off.

f) Curing environment

The components to be cured can be

cured under different environments. The environment to be set depends on the application and can be found in table 1.

g) Open system

In the PCU LED, 10 freely programmable storage spaces can be used.

Table 1: Curing environment for different uses

Environment	Use	PCU LED vacuum	PCU LED nitrogen
Vacuum	Curing of laboratory products	✓	✓
Nitrogen	Curing of medical devices		✓

2. Technical data

Weight of the overall package	11.1 kg
Weight of the device	9.3 kg
Dimensions of the overall package	H 220 x W 455 x D 370 mm
Dimensions of the device	H 110 x W 389 x D 276 mm
Polymerization chamber	H 65 x W 150 x D 150 mm
Voltage	100–240 V
Nominal frequency	50–60 Hz

3. Declaration of conformity

According to Machinery Directive (2006/42/EG), Low Voltage Directive (2014/35/EU) and EMV Directive (2014/30/EU)

Herewith we declare that the machine described below in its conception and design and in the shape delivered by us is in accordance with the fundamental requirements for safety and health as prescribed in the applicable EC directives. In the case of any change or modification of the machine, not authorised by us, this declaration becomes invalid.

Denomination of the machine: PCU LED/PCU LED N₂
^{REF} 4316/4317

Type of machine: Light polymerization unit

Applicable EC directives: 2006/42/EC According to Machinery Directive
 2014/35/EC Low Voltage Directive
 2014/30/EC EMV Directive

Applied harmonized standards:

2006/42/EC According to Machinery Directive

EN ISO 12100:2010 Safety of machinery – General principles for
 + Corrigendum to design – Risk assessment and risk reduction
 DIN EN ISO 12100:2011-03 (ISO 12100:2010); German Version EN ISO
 12100:2010 + Corrigendum to DIN EN ISO
 12100:2011-03

2014/35/EC Low Voltage Directive

EN 60204-1:2006/A1:2009 Safety of machinery – Electrical equipment
 + Corrigendum to EN 60204-1:2006 of machines – Part 1: General requirements
 (IEC 60204-1:2005/A1:2008); German Version
 EN 60204-1:2006/A1:2009 + Corrigendum to
 EN 60204-1:2006; German Version CENELEC-
 Cor.:2010 to EN 60204-1:2006

EN 61010-1:2010 Safety requirements for electrical equipment
 for measurement, control and laboratory use
 – Part 1: General requirements (IEC 61010-
 1:2010 + Cor.:2011); German version EN
 61010-1:2010

2014/30/EC EMV Directive

EN 61000-3-2:2014 Electromagnetic compatibility (EMC) – Part 3-2:
 Limits – Limits for harmonic current emissions
 (equipment Input current ≤ 16 A per phase)
 (IEC 61000-3-2:2014); German version EN
 61000-3-2:2014

EN 61000-3-3:2013

Electromagnetic compatibility (EMC) – Part
 3-3: Limits – Limitation of voltage changes,
 voltage fluctuations and flicker in public low-
 voltage supply systems, for equipment with
 rated current ≤ 16 A per phase and not sub-
 ject to conditional connection (IEC 61000-3-
 3:2013); German version EN 61000-3-3:2013
 Electromagnetic compatibility (EMC) – Part
 3-3: Limits – Limitation of voltage changes, vol-
 tage fluctuations and flicker in public low-vol-
 tage supply systems, for equipment with rated
 current ≤ 16 A per phase and not subject to
 conditional connection (IEC 61000-3-3:2013);
 German version EN 61000-3-3:2013

EN 61326-1:2013

Dreve Company Group
 Max-Planck-Straße 31 · 59423 Unna/Germany
 Tel.: +49 2303 8807-0
 Fax: +49 2303 82909
 E-Mail: info@dreve.de
 www.dreve.com

4. Safety instructions

Attention! Please read these notes carefully before the connection and commissioning of the device. The operational safety and functioning of the device can only be ensured if both the general safety and accident prevention regulations as specified by the legislator and the safety instructions in the operating instructions are observed.

1. The device must only be used in accordance with the present operating instructions. We do not assume any liability for damages which result from improper use and/or faulty operation.
2. Only operate the device at a socket with a protective conductor. Never touch the plug with wet hands.
3. The device must be placed on a stable and level surface which is sufficiently

stable and loadable for the operating weight of approximately 9.5 kg.

4. Do not insert any objects into the device apart from the polymerization chamber.

5. Only accessories and spare parts may be used which are approved by the manufacturer.

6. Do not lock the operating and safety elements of the device.

7. Please check the device to ensure its proper condition and operating safety before each operation. If its condition is not flawless, the device must not be used and must be labelled accordingly.
8. Before the cleaning and maintenance of the device or before replacing parts, the power plug must be pulled out.
9. Signs and stickers must always be kept in an easy-to-read condition and must not be removed.
10. The opening of the device and repairs may only be performed by authorized professionals.
11. Only use nitrogen as protective gas. Do not use flammable or toxic gases under any circumstances. Please observe the general safety and accident prevention regulations when handling nitrogen.
12. Only use approved gas pressure bottles.
13. When connecting the gas pressure bottle to the device, be sure to check the tightness of the feed tube. Please make sure that the gas pressure bottle cannot tip over.
14. Never connect the gas bottle without using a pressure reducer; the device may be operated with a maximum pressure of 5 bar (70 psi).
15. Please close the gas supply, when the device is not used. Do not operate the device unattended.
16. Be sure to provide for sufficient ventilation in the work area.
17. The pressure of the gas line must not exceed 5 bar (70 psi).
18. Do not keep any highly inflammable substances within the immediate vicinity of the device.
19. Authorized operators: The operator of the device must make the operating instructions accessible to the operator and must also make sure that s/he has read and understood them. Only then may the operator put the device into operation.
20. To prevent the ingress of water into the device (e. g. splash water), the device should be installed in a dry environment.
21. Please disconnect the device from the mains when it is not used.
22. To avoid the risk of getting caught, only touch the lid at the handle.
23. Do not operate the device if the housing of the device is not closed and screwed on.
24. Unauthorized modifications and changes are inadmissible for safety reasons.
25. The prescribed operating and maintenance conditions specified in these operating instructions must be complied with. Please observe the general accident prevention regulations when working with the PCU LED.

5. Commissioning

Please ensure before the commissioning that no damage has occurred during transport. Check the housing for external damages and control the glass plates on the floor and at the lid of the polymerization chamber. The lid is secured for transport reasons. The transport locks must be removed before the commissioning of the device. Please ensure that the indicated voltage on the type plate corresponds to the voltage of the power source. Insert the power cable into the appliance socket and create a network connection with the connection cable.

Attach the protective gas supply to the connection of the pressure reducer of the gas pressure bottle with the use of

the hose clamp. Establish a connection to the device with the quick-lock coupling. The pressure must not exceed 5 bars (70 psi).

Scope of delivery:

The scope of delivery varies depending on the equipment:

1x PCU LED

1x Accessory package

The accessories consist of:

3x Power cable

1x Operating instructions

1x Gas connection line kit including quick coupling and 2 hose clamps (only for nitrogen option)

1x Post curing insert

(only for nitrogen option)


6. Information on the use of the protective gas (PCU LED (N₂))

The use of the protective gas allows for the completion of the polymerization in an atmosphere without oxygen. The protective gas replaces the air in the polymerization chamber and thus reduces the oxygen content in the polymerization space of the PCU LED, which leads to the hardening of components without an inhibition layer. The protective gas supply is automatically controlled when the polymerization process starts. The flow rate is pre-set; the adjustment can be changed at the control needle valve at the back of the device (please observe the safety instructions).


The outlet pressure of the bottle can be mounted to the pressure reducer of the gas pressure bottle so that the flow is changed without changing the control needle valve. The gas flow is factory-set. The use of a flow meter is recommended for monitoring the gas flow during operation. (Please refer to the instructions for use of the relevant material for exact information on type and flow rate.) The connection pressure must be set to 5 bars.

7. Functional elements

Overview

Rear view


8.1. Operation: Starting the device

The software of the PCU LED is subdivided into two levels: The first level is the user level, where the operator can select and confirm different options. The second level is the process owner level, where the process owner can perform different device configurations.

- (1) Switching on the power plug at the back of the device brings you directly to the user level.
- (2) Switching on the power switch at the back of the device and the simultaneous pressing of the operating button brings you to the process level.

8.2. Functions of the rotary button

Selection

Confirmation

Selection of the sub-menu

8.3. Operation: User level

Start screen after switching on the device


Welcome screen

Firmware version

Selection of the material parameters


If a parameter set is selected, it will be highlighted on the display. Additional confirmation brings you to the selected parameter-set in the status window.

Status window


Additional confirmation starts the exposure process.

- 1 – Programme name
- 2 – Return button
- 3 – Exposure time
- 4 – Pressure (mbar)
- 5 – Exposure icon (ON/OFF)
- 6 – Exposure atmosphere (air/vacuum/nitrogen)

Ventilation of the exposure room


Wait until the vacuum is released so that e. g. the sealing for the exposure chamber is not destroyed!

Final process screen

After the exposure process, the result is shown together with the used parameter sets.

8.4. Operation: Process level


Access to the configuration level

Switching on the power switch at the rear of the device and simultaneously pressing the control button takes you to the process level. (See chapter 8.1)

Entry of the access code 217 (Status 03 / 2017)


Please confirm the code after entry. The access code cannot be freely selected and can also change in the case of a firmware update.

Configuration level in overview

1

1 – Orientation line:
A vertically oriented narrow bar at the left and the right margin of the selected window is provided for guidance purposes


Selection of the process parameters

10 programme blocks available for programming

Green – activated for the user level


White – not activated for the user level

Setting the process parameters – press the rotary button > 5 seconds

Setting the process parameters


- 1 – Return button
- 2 – Programme name (freely definable by turning the rotary button)
- 3 – Intensity setting for LED sources (freely definable by turning the rotary button)
- 4 – Exposure medium
- 4.1 Green – nitrogen active
- 4.2 Green – vacuum active
- 4.1 and 4.2 not marked – air active
- 5 – Exposure time

Setting the LED intensities


The intensity can be shown in percent per panel. There are 3 LED spots on each panel.

Attention: If the predetermined values are modified for the appropriate material, the final material properties indicated by the manufacturer can be changed.


Checking the vacuum function

The vacuum function is checked on the basis of a special vacuum value being reached after a defined time. The result is shown on this screen after the test.

Time setting

Setting the date and time for the time stamp of the process recording which is shown, inter alia, in the "final screen" in the user level.

Factory setting

When reverting back to the factory settings, the individual parameter settings are deleted and the device is reset to the delivery state.

8.5. Placing the components

With the nitrogen version, the supplied post-curing insert is placed in the exposure chamber and the components are

positioned on the glass plate of the post-curing insert.

9. Maintenance

Maintenance work to be performed by the user

Maintenance unit	Maintenance frequency
Clean the glass panes in the bottom and the lid of the polymerization chamber with isopropyl alcohol	* daily
Before commissioning check whether all connections of the device are positioned correctly: <ul style="list-style-type: none"> • gas connection • mains connection 	* daily
Cleaning of the device with a mild detergent	* weekly

Maintenance work to be performed by the service commissioner

Maintenance unit	Maintenance frequency
Check whether the LED panels must be or changed	* weekly with program LED-Test
Check whether the LED panels must be calibrated	* monthly
Check whether the gas valve works reasonably	* weekly
Check the lid closing function	* weekly
Check whether all seals are still functional and/or fit correctly	* monthly

* Recommendation, actual need depends on use intensity and other requirements.

10. Monitoring

The different process parameters are recorded by means of the sensors used in the PCU LED. These are shown in the final screen in the user level.

Explanation of the parametric representation in the final screen in the user level

Presentation	Meaning
216	Sequence number (number of exposure jobs)
01/01/2000	Date
01:17:35	Time
040	LED – intensity of the upper panel in %
040	LED – intensity of the lower panel in %
0800	Exposure time [place 1 and 2: minutes; place 3 and 4: seconds]
1	Vacuum ON [1 – ON; 0 – OFF]
0	Nitrogen OFF [1 – ON; 0 – OFF]

Breakdown of the error codes in the final screen in the user level

Code	Status	Background display
00	Accurate work flow	Green
01	Pump runs during the exposure process	Red
02	Pump cannot draw any vacuum	Red
03	Lid opened	Red
04	No nitrogen flow	Red
05	Run time of program = 0	Red
06	LED error in program LEDTest	Red

11. Troubleshooting plan


Code	Defect	Cause	Troubleshooting
00	Accurate work flow	Not applicable	Not applicable
01	Pump is running during the exposure process	Cover gasket of exposure chamber leaky	Replacement of the lid sealing
		Hoses of the vacuum system leaky	Replacement of the hoses

Code	Defect	Cause	Troubleshooting
02	Pump cannot draw any vacuum Cutout after 150 seconds	Cover gasket for exposure chamber leaky	Replacement of the lid sealing
		Hoses of the vacuum system leaky	Replacement of the hoses
03	Lid is opened LED do not work	Vacuum pump defective	Exchange of the vacuum pump
		Lid not closed	Close lid
04	No nitrogen flow	Lid control (micro switch) defective	Exchange lid control (micro switch)
		LED defective	Exchange of the LED
05	run time of program=0	Nitrogen valve at the gas cylinder not opened. Warning: The device is only o be operated with a maximum pressure of 5 bar (70 psi)	Open nitrogen valve gas cylinder not opened.
06	LED error in program LEDTest	no run time of program was set	set a run time > 0
06	LED error in program LEDTest	LED defective	LED has to be changed by trained personnel

Further

Device does not show any function	Network connection not correct	Control network connection
	Main switch off	Switch on main switch
	Main switch defective	Exchange main switch

12. Exploded view

13. Spare parts list

	REF
1 button cell	51090CR2032
2 PCU LED control	521022
3 power supply	521023
4 glass plate	501024
5 mirror plate, rear side lower part	521025
6 mirror plate, rear side cover	521026
7 mirror plate, sideways right	521027
8 mirror plate, sideways left	521028
9 rotary button	521029
10 LED panel	521030
11 mirror plate, front side	521031
12 adaptor for rotary button	521044
13 vacuum seal	521050
14 voiframe - display cap	521057
15 cap for measuring feed-through	521063
16 vacuum pump	521064KM32
17 UK power supply angled	51284
18 plate	521065
19 cable harness PCU LED with USB bushing	521066
20 angle piece	521075
21 compression spring	521078
22 cold device power cable	521512
23 locking screw	521083
24 magnetic valve	521560
25 gas compression spring	521565
26 encoder	521566
27 power input + filter	521568
28 fabric tube	2004
29 microfuse	50908
30 choke valve	50937
31 sealing ring	51012
32 mecanyl pipe	51030
33 threaded connection	51050
34 enclosure stand	51079
35 USA power supply straight	51283
36 Glass plate for PCU LED insert	521092
37 Frame of PCU LED insert	521091

14. General information

Repair and maintenance work at the PCU LED may only be performed by specialists. The power supply to the device must be disconnected in any case.

The company Dreve will be responsible for the safety, reliability and performance of the device, if

- the device is only used for the appropriate application purpose.
- the operation is only carried out by professionally trained personnel.

- the device is used in accordance with these operating instructions.

- extensions, new settings, changes or repairs are only performed by authorized persons.

Since we at Dreve are constantly working on the further development of our products, we cannot exclude the possibility of technical changes.

15. Warranty

The manufacturer warranty is 2 years from date of delivery. The warranty right (liability of defects) is not limited by the warranty and is regulated in the general terms and conditions of the company Dreve. The guarantee includes the assumption of repair and shipping costs in the case of material and manufacturing defects. The warranty claims do not refer to natural wear nor to damages which are the result of the transfer of risk due to faulty or unintended use, excessive loading, unsuitable equipment and any such chemical, electrochemical or electrical influences which are not provided in accordance with the contract. Warranty claims are excluded in the case of changes or

repair and maintenance work which have been incorrectly completed by the purchaser or third parties as well as in the case of the breaching of the seals.

All technical changes to the machine are prohibited and will result in the loss of the warranty claim. Repairs may only be performed by trained personnel of the company Dreve. For cases of repair, please contact the company Dreve, Max-Planck-Str. 31, 59423 Unna / Germany; Phone: +49 2303 8807-0.

After the sending-in of the product and the submission of the sales receipt, the product will be repaired free of charge or replaced by a comparable product.

16. Decommissioning of the product

The device is subject to the provisions of the Elektro- und Elektronikgerätegesetz [German Electrical and Electronic Equipment Act], the return of this device is guaranteed by the public bodies. This regulation applies to all member states of the European Union.

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