When embedding FotoDent® cast, please note the manufacturer's data of the used embedding materials. To ensure that the resin combusts without any residue, a holding time of 60 minutes at 800–900 °C is recommended.

To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and colour variations of the FotoDent® cast material.

Possible dangers

Please note the safety instructions for processing FotoDent® cast.

Hazard statements: H315 Causes skin irritation. H319 Causes serious eye irritation. H317 May cause an allergic skin reaction. H335 May cause respiratory irritation.

Precautionary statements: P261 Avoid breathing dust/fume/gas/mist/vapours/spray. P280 Wear protective gloves/protective clothing/eye protection/face protection. P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P405 Store locked up. P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains: Diurethandimethacrylate, triethylene glycol dimethacrylate, phenyl bis(2,4,6-trimethylbenzoyl)-phosphineoxide

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® cast packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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Intended use

FotoDent® cast is a light-curing resin for the production of cauterization profiles for casting technology using 385 nm/405 nm LED based stereolithography. For post curing of FotoDent® cast built parts, the light curing units FotoDent® flash and PCU LED are recommended, which are optimal balanced to the production process.

Processing

- During the application of FotoDent® cast we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- FotoDent® cast should be well shaken approx. 1 hour before use.
- Carefully pour FotoDent® cast into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.
- Make sure that FotoDent® cast is tempered in your production unit up to approx. 23 °C.
- Select the product specific parameter settings for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 10 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended.
- Post-curing:

PCU LED	6 minutes, 80 % light power (vacuum)
FotoDent® flash	8 minutes + 1 minute gas pre-flow
- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® cast. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1 hour with FotoDent® cast.