- For post-curing of the parts, please select the parameter settings for your process. You can find more qualified printers on our website www.dentamid.dreve.de.

Post curing	
PCU LED N ₂ Rapid Shape D, DII / Asiga MAX UV Sonic XL 4K (qualified by Dreve) / Sonic XL 4K 2022 (qualified by Dreve)	6 minutes, 80 % light power (vacuum)
FotoDent® flash Rapid Shape D, DII Asiga MAX UV	8 minutes + 1 minute gas pre-flow 10 minutes + 1 minute gas pre-flow
RS Cure	According to respective parameter set for FotoDent® cast.

- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® cast. It is recommended to regularly take out the container of the production unit, to homogenise and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hour with FotoDent® cast.

When embedding FotoDent® cast, please note the manufacturer´s data of the used embedding materials. To ensure that the resin combusts without any residue, a holding time of 60 minutes at 800–900°C is recommended.

To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and/or colour variations of the FotoDent® cast material.

Hazards identification

Hazard and safety notes can be found on the product label and the corresponding safety data sheet.

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® cast packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

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Product description

Light-curing material for the production of residue-free burn-out moldings for dental casting technology.

Processing

- During the application of FotoDent® cast we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.
- FotoDent® cast should be well shaken approx. 1 hour before use.
- Please see further steps within printer manual.
- Carefully pour FotoDent® cast into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.

Processing temperature	
Sonic XL 4K (qualified by Dreve) / Sonic XL 4K 2022 (qualified by Dreve)	23 ± 3 °C
Rapid Shape D and DII printer	23 ± 3 °C
Asiga MAX UV	35 ± 3 °C

- Select the product specific parameter settings for FotoDent® cast for your process.
- When the building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 2x6 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended. Alternatively, cleaning with isopropanol can be carried out in a separate beaker under pouring conditions for approx. 2x6 minutes. Recommendation: Beaker and magnetic stirrer, laboratory shaker or ultrasonic cleaner. In case of insufficient cleaning (such as shiny areas on the object after drying) repeat cleaning process.