## FotoDent® guide

(ES)

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## FotoDent® guide

### Intended use

Material on methacrylate resin basis for DLP systems with 385nm resp. 405 nm-LED for manufacturing of dental surgical drill guides for tooth-supported and/or teetheridge-supported and surgical invasive application.

**Indications for use:** Manufacturing dental drill guides for placing dental implants.

**Indication:** Generative manufacturing of drill guides using 3D printing. This serves the dentist as a guide rail during drilling.

**Contraindication:** If an acrylic allergy exists, an allergic reaction may occur.

**Patient target group:** Patients with a single-tooth gap, partially and/or completely edentulous upper and/or lower

**Intended users:** The dental technician or specialist staff in the dental laboratory is responsible for the production of the drill guide in 3D printing. The dentist uses the finished drill guide on the patient.

### 1.Requirements/equipment

**Printer:** Rapid Shape: D and DII series / Asiga: MAX UV; PRO 4K / Dreve: Sonic XL 4K (qualified by Dreve) / Sonic XL 4K 2022 (qualified by Dreve)

**Software:** Autodesk Netfabb®, Composer, Dreve Elements

**Cleaning:** Isopropanol, ultrasonic cleaner

**Post curing unit:** PCU LED N<sub>2</sub>, FotoDent® Flash, RS cure.

Please contact us if you are using any other type of equipment.

### Basic Material Properties

Please refer to Technical Data Sheet (TDS). Available on request.

### 2. Process description

#### 2.1 Preparation

• Alignment of drill guide: vertical, supports, connection of supports to a baseplate.

• During processing of the liquid material, we recommend wearing personal protective equipment e. g. suitable gloves, safety goggles etc.

• Carefully pour the material into the pre-determined container of the production unit.

• Please see further steps within printer manual.

• Remove all bubbles with a cleaned object.

#### 2.2 Printing process

• Processing temperature: See overview.

• Select the corresponding build style/INI file for FotoDent® guide in the respective color in the printer software (material library). Make sure the software used is up to date.

Rapid Shape DII: FotoDent guide 385nm #D35650
Rapid Shape D: FotoDent guide 405nm #D35600
Asiga: FotoDent guide\_D35650.ini
Dreve: Dreve FotoDent guide405

• Start the printing process.

#### 2.3 Processing after the printing process

• It is recommended that the printing result be further processed directly after printing (see following steps).

• After the platform is taken up a dripping off time of approx. 10 minutes isrecommended.

#### 2.4 Cleaning

• Cleaning and post curing incl. supports and baseplate to avoid warpage.

• Clean the drill guide in isopropanol medium (purity ≥ 97 %) in the ultrasonic bath. Fill two suitable containers with isopropanol.

1. Pre-clean with isopropanol for 6 minutes in the first suitable container.

2.Blow out the drill guide using compressed air, pay attention to any additional holes.

3. Perform the main cleaning step with isopropanol for 6 minutes in the second suitable container.

4. Blow out the drill guide using compressed air, pay attention to any additional holes.

#### 2.5 Post-curing

See overview post-curing.

• After post-curing, the supporting structures and base plate must be removed.

### 2.6 Finishing

• The surface can be smoothed by sanding or polishing.

• Fitting of suitable drill sleeves (from other manufacturers) into the drillguide. Please refer to the instruction for use of the supplier/manufacturer for the adjustment of the drill sleeves.

### 3. Important notes

• Select FotoDent® guide material that matches the wavelength of the used printer (385nm or 405 nm).

• The specifications and biological safety have been qualified using the aforementioned printer, the associated software and the process parameters indicated.

• To avoid deterioration of the material quality, protect the liquid material from exposure to light. Variations from the described manufacturing process may impair biological safety, lead to modified mechanical characteristics and/or color variations of the material.

• Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. However, thanks to the low viscosity, it is possible to filter FotoDent® guide. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hours with FotoDent® guide.

### 4. Cleaning and disinfection

For cleaning and disinfection, we recommend the use of ethanol (purity ≥ 60 %) or a 0.2 % chlorhexidine solution, each with a contact time of 5 minutes, then allow to dry well.

### 5. Sterilisation

Do not use thermal methods for sterilisation, otherwise deformation of the surgical guide cannot be ruled out. FotoDent® guide cannot be sterilised.

#### 6. Serious incidents

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

### 7. Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® guide packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

#### 8. Hazards identification

Hazard and safety information can be found on the product label and/or in the corresponding safety data sheet.

#### 9. Disposal

Disposal of contents/container in accordance with local/regional/national/international regulations and according to the safety data sheet.

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**Drucker:** Rapid Shape: séries D and DII / Asiga: MAX UV; PRO 4K / Dreve: Sonic XL 4K (qualified by Dreve) / Sonic XL 4K 2022 (qualified by Dreve)

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