

## Technical Data Sheet

**Product description:** Cleaning tablet for external cleaning of dental splints

**Manufacturer:** Dreve Dentamid GmbH  
Max-Planck-Straße 31  
59423 Unna / Germany

### Composition:

Pos	Chemical term	CAS-No.
1	Potassiumperoxomonosulphate	70693-62-8
2	Sodiumbicarbonate	144-55-8
3	Anionic Tenside	151-21-3
4	Sodiumpercarbonate	12630-89-4
5	Citric acid	77-92-9
6	Food couring	
7	Mint oil	

### Properties:

Diameter: 20 mm  
Thickness: 5 mm  
Weight: 3,0g ± 0,15g  
Colour: tinged blue white  
Odour: peppermint-like  
Period of decomposition: < 5 min  
ph-Value: 7,0 – 9,0  
Packaging: Alu-Blisterfolie

Dreve Dentamid GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 8807-55

[info@dreve.de](mailto:info@dreve.de)

[www.dreve.com](http://www.dreve.com)