

# Sonopal<sup>®</sup> soft Lack E

## Technical Data Sheet

**Innovation  
MediTech GmbH**

Max-Planck-Straße 31  
59423 Unna / Germany  
Tel.: +49 2303 8807-42  
Fax.: +49 2303 82909

Product description: Light curing lacquer for BTE and ITE hearing systems.

Manufacturer: Innovation MediTech GmbH  
Max-Planck-Straße 31  
59423 Unna / Germany

Composition:

| Pos | chemical term                |
|-----|------------------------------|
| 1   | initiator system             |
| 2   | methacrylate                 |
| 3   | acrylates                    |
| 4   | aliphatic urethane acrylates |

Form: fluid

Density: 1.00 – 1.02 g/cm<sup>3</sup> at 20 °C

Viscosity: 0.02 – 0.05 Pa s at 23 °C

These data result from measurements which are executed from time to time in view of our quality management system.  
This document is valid without any signature.

Rev. 4 / 10.2024  
Page 1 of 1  
TDS\_Sonopal soft  
Lack E\_GB